



CE

Tabele of contents

1. Introduction	3
2. Description of product	3
3. Application	4
4. Contradictions.....	4
5. Usage	4
6. Safety of use - WARNING !	5
7.1. Remote control	8
7.2. Headrest adjutment.....	8
7.3. Knee protectors.....	9
7.4. Footrest adjustment	10
7.6. Backrest.....	12
7.7. Wheels.....	12
7.8. Replacement tires	13
7.9. Batery charging.	14
7.10. Problem solving.....	15
7.11. Batteries.....	15
7.12. Transport	16
7.13. Maintenance	16
8. Technical paramets	17
9. Symbols	18
10. Contact details	19

**DO NOT USE THE WHEELCHAIR FOR THE FIRST TIME WITHOUT
COMPLETETY READING THE USER'S MANUAL.**

1. Introduction

Your Baby Hero wheelchair combines the most advanced state-of-the-art components with modern, attractive styling. We are sure that the design features and trouble-free operation will add convenience to your daily living and ensure complete satisfaction.

Your safety is most important to us. **Please read and follow all the instructions outlined in this manual before you attempt to operate your wheelchair for the first time.** These instructions were produced for your benefit. Your understanding of these instructions is essential for safe operation of your new Hero wheelchair. Manufacturer is not liable for damage to property or personal injury arising out of the unsafe use of wheelchair. Manufacturer is also not liable for any property damage or personal injury arising out of the failure of any person and/or user to follow the instructions and recommendations set forth in this manual or any other instructions or recommendations contained in other Hero related literature issued by manufacturer.

This manual presents the newest current technology standards and product information. As we constantly improve technology and our products, you may find some differences between explanations and illustrations presented in this manual and product you have purchased. However, it does not affect regular use of your wheelchair. If you have any questions which you can't answer while using the wheelchair (safely, or while you think according to the instruction, and/or when the recommendation of the manual cannot be carried out safely) please, contact your local dealer.

Please, check and read the warranty card carefully. Please, make sure that it has been properly filled in and stamped. Keep it in case you will have to use it in the future. The warranty card is an integral part of the product you purchased, if you sell your wheelchair, do not forget to pass it (together with user's manual) to the next owner. Always have the user's manual when you use your wheelchair in case of any problems. If you are not able to solve the problem or if you are not sure if you can use the wheelchair safely, please, consult your dealer.

2. Description of product



Illustration1.

1. Headrest
2. Handgrips
3. Backrest
4. Armrests
5. Guard rail
6. Positioning belts
7. Seat
8. Knee protectors
9. Heel loop
10. Footplate
11. Stabilizing footer
12. Caster front wheels
13. Rear wheels



Illustration 2.



Illustration 3.

3. Application

Baby Hero Manual standing wheelchair can be used in patients with congenital or acquired paralysis of the lower limbs, muscular dystrophy, stroke and amyotrophic lateral sclerosis. It can be used also for the prevention of pressure ulcers.

4. Contradictions

Contraindications to a sitting position and upright. Patients with mental disorders should not remain in the standing wheelchair unattended. During lifting phase observe the reactions of the patient. If you notice any adverse reactions, eg. turn pale, sweating, fainting should immediately stop verticalization and return to the starting position.

5. Usage

Baby Hero Manual standing wheelchair thanks to standing option:

- Prevents joint contracture of lower limbs;
- Improves bowel & bladder's functions;
- Prevents skin abrasions and bedsores arising from long sitting position;
- Prevents osteoporosis – improves calcification, especially of lower limbs;

- Improves self-confidence and reduces the fear of assuming standing position;
- provides body weight support training, prepares the patient to assume standing position by himself;
- improves cardiopulmonary functions;
- allows to perform every-day activities in standing position which increase sense of personal independence;
- Increase adaptation ability.

6. Safety of use - WARNING !



- Not following the instructions of the user's manual may result in personal injury or damage to the product.
- Failure to follow the safety instructions in the manual may result in injury.
- To increase the comfort and safety during using wheelchair, should follow suggestions given in the instructions.
- During transport the power cable connected to the battery should be disconnected to protect the life of the battery. Before using the wheelchair, remove the white cap and plug the battery.



Illustration 4.

- Do not modify your wheelchair in any way without manufacturer's authorization. Unauthorized modifications may result in personal injury and/or damage to the product.**
- Be particularly careful when mounting, dismounting, verticalization and folding the lack of attention may result in finger trap by moving parts.**
- Before the first use of the standing frame is necessary to carefully read the entire instructions for use.**
- When using the wheelchair for the first time, assistance of other person is necessary.
- As you begin using your wheelchair during daily activities, you will probably encounter situations in which you will need some practice. In such situations, take your time and you will soon be in full and confident control as you manoeuvre through doorways, on and off elevators, up and down ramps, and over moderate terrain.
- Transferring onto and off your wheelchair requires good sense of balance. If you start learning how to get on/off the wheelchair, assistance and help of other person is necessary.
- If you need to use the positioning belts, make sure they have been property fastened and adjusted to ensure maximum safety
- Please, remember that Hero should not be used on public street and road – it can be difficult for traffic to see you when you are seated on your wheelchair.
- Never driver down the stairs or escalators.
- Be careful while driving over doorsteps or other obstacles.

- o) Do not bend, lean or reach for objects while sitting on the wheelchair – the wheelchair may tip which can cause product damage and personal injury.
- p) Always protect batteries from freezing. Never charge frozen batteries – it may lead to battery damage and personal injury.
- q) If you take prescribed or over-the-counter medication or if you have certain physical limitations, please,
- r) consult your physician. Some medications and limitations may impair your ability to operate your wheelchair safely.
- s) Never operate your scooter while you are under influence of alcohol, as this may impair your ability to
- t) drive safely.
- u) Do not put your overall body weight to the footrests – the wheelchair may tip and cause personal injury.
- v) Hero wheelchair is designed for one person only, you must not carry any passengers!
- w) Nie należy pozostawiać Pionizatora wraz z użytkownikiem w bezpośrednim pobliżu otwartego źródła ognia.
- x) Do not operate your wheelchair when the temperature falls below 0°C, in icy or slippery conditions or on salted surfaces (i.e., walks or roads). Such use may result in an accident, personal injury, or adversely affect the performance and safety of your scooter.
- y) Do not operate or store your wheelchair in a place where it may be exposed to inclement weather condition such as rain, snow, mist, and below freezing temperatures. Using the wheelchair in such conditions can damage the electronics and potentially result in loss of control. Before using your wheelchair, make sure that the water has dried up.
- z) Be aware that the surface temperature of the upholstery and the frame parts may be increased after exposure to the sun or other heat source. It may cause burns parts of the body.
- aa) Do not move wheelchair by stairs. If you have to overcome the stairs you have to be transported in a different way. Wheelchair must be transported separately.
- bb) Do not attempt to lift the standing frame holding it by moving parts or possible to remove, because it can cause accidents and injuries.
- cc) Never use the wheelchair as a seat for the transportation of people in cars, buses or any motor vehicles.
- dd) The wheelchair is designed for people up to 100 kg
- ee) Wheelchair must be easy to fold and to unfold. In case of problems with folding, you should re-read the instruction or contact the manufacturer.
- ff) Brakes cannot be slow down the wheelchair during ride, You should use them only to prevent the unintended movements
- gg) Both brakes should be properly adjusted. Incorrect setting one of the brakes, results with blocking only one rear wheel. On large slopes, it can cause difficulty to control the rotation the wheelchair around the wheel stopped.
- hh) Regularly check the tire pressure. With too low tire pressure than recommended, the brake performance can be significantly reduced!
- ii) All activities, such as getting in, getting off, changing the position of the patient, stop wheelchair require block the brakes.
- jj) Do not leave the patient in the wheelchair (even with turned on both brakes) on the descents, ramps, landings, etc.
- kk) It is necessary to check after each wheel mount if wheel axle connector lock works properly. While mounting should hear “click” of the latch. Pull strongly the wheel to check if it is installed correctly.

- ll) Using the parking brake while driving may result with an accident with all its consequences.
- mm) After each wheelchair adjustment, it is always necessary to check the correct operation of the brake.
- nn) Oil, grease, water or other chemicals, may cause that brakes may have reduced effectiveness.
- oo) Strength of the breaking intensity may be reduced after long exploitation of the wheelchair.
You are allowed to adjust breaks yourself. In order to do that you must loosen the screws locking the position of the break, move the mechanism of the break towards the wheel so that the wheel will not rotate when break will be on again (try the adjustment with the user sitting on the wheelchair). After final adjustment secure back the mechanism with the screws.
- pp) It is extremely important to keep all wheelchairs components clean in order to maintain the product parameters.
- qq) It is unacceptable to use the wheelchair after discovering any visible or detectable damages or when the wear of some parts may cause an accident. All the repairs need to be performed by a skilled personnel and during the warranty period only by the manufacturer's service.
- rr) In case of wheelchair fitted with pumped wheels You should not exceed the maximum permissible inflation pressure of the wheels, the value is indicated on the tire wheels.
- ss) During use, handing of the product and during its assembly, regulating mechanisms, it might happen trapping and / or squeeze user's or caretaker's body part in holes / gaps between the elements. Perform these steps with extreme caution. After adjustment, stabilize the position by carefully tightening the nuts / bolts.

RAMPS , OBSTACLES

- Approach the ramp/obstacle slowly.
- Do not drive up/down the incline of angle more than 10°.
- Avoid driving up/down long inclines.
- Do not leave your wheelchair on any incline – no matter whether you seat or stand next to it.

DRIVING OUTDOOR

Baby Hero and Hero wheelchairs have been designed to provide you with maximum comfort of use and safety when driving in standard conditions – on flat and dry surface.

- Avoid high grass – it may become tangled in the running gear.
- Avoid loosely packed gravel and sand.
- If you feel unsure about a driving surface, avoid that surface.

Please, follow the instruction below if you need to transport the wheelchair between the floors and the elevator is not available. The user needs to get off the wheelchair.

- Lift the wheelchair (attention! Do not hold its removable parts!)- remember to keep your back straight and to bend your knees.
- Transfer the Hero up/down the stairs.
- Put the wheelchair only on the flat surface.

Safety belts

Always remember to fasten the thigh and thoracic safety belt immediately after taking a seat in the wheelchair.



Illustration 5.





Illustration 6.

7. Everyday use and maintenance

7.1. Remote control



Illustration 7.

1. OFF/ON LED light.
2. Press and hold the button  to assume standing position.
3. Press and hold the button  to return to sitting position.
4. Battery condition meter .

The battery condition meter indicates the battery status. There are three LED lights:

- Green means fully charged batteries
 - Red means that batteries are charged only in 30 % and that you need to charge them IMMEDIATELY
- Red LED light also flashes when the batteries are being charged.

7.2. Headrest adjustment

- In order to adjust the headrest turn the rear button in anticlockwise direction. Choose desired height and turn the rear button in clockwise direction. (Illustration 8.)




Illustration 8.



Illustration 9.

- In order to adjust vertical inclination of headrest, release the screw (see the picture) using hexagon spanner. Choose the desired angle and fasten the screw tightly. (Illustration 9.)
- Release the butterfly screw of the headrest – choose the desired angle and fasten the butterfly screw. (Illustration 10., 11., 12.)



Illustration 10.



Illustration 11.



Illustration 12.

7.3. Knee protectors

- Unlock the blockade according to the picture – move the footrest outside. (Illustration 13.)



Illustration 13.



Illustration 14.



Illustration 15.

- If you unlock the blockades on both sides, you will be able to remove the knee protectors. (Illustration 14.)
- Adjust length of heel loop. (Illustration 15.)
- In order to adjust height and depth of knee protectors, release two screws placed on both sides of the wheelchair. (Illustration 16.)




Illustration 16.

7.4. Footrest adjustment

- Unfasten the heel loop. (Illustration 17.)
- Lift the footrests. (Illustration 18.)



Illustration 17.



Illustration 18.

- Unlock the footrests according to the picture. (Illustration 19.)



Illustration 19.



Illustration 20

- Move the footrests outside. (Illustration 21.)
- Lift and remove the footrests.
- In order to adjust footrest height, release the screws located at the back of the footrests. (Illustration 22.)



Illustration 21.



Illustration 22.

- There can be lumbar cushion under the backrest (optional). (Illustration 23., 24.)



Illustration 23.



Illustration 24.

7.5. Armrests and guard rail

There are screws and additional holes which allow you to adjust armrests height. Release the armrests, choose desired height and fasten the screws tightly. (Illustration 25.)

- In order to lift the armrests, release the screws, choose desired angle and fasten the screws. Always remember to choose the same height for both armrests. (Illustration 26.)



Illustration 25.



Illustration 26.

- Additional guard rail can be installed. (Illustration 27.)
- In order to remove the guard rail, push the switch placed on the rear part of the armrest and remove the guard rail. (Illustration 28., 29.)



Illustration 27.




Illustration 28.




Illustration 29.

- In order to facilitate getting on and off the wheelchair, the guard rail can be moved aside. (Illustration 30., 31.)



Illustration 30.



Illustration 31.

7.6. Backrest

- In order to remove the backrest, release the screws placed at the back of the backrest and pull up the backrest. (Illustration 32., 33.)



Illustration 33.



Illustration 33.

7.7. Wheels

- In order to block the wheels, push the hand brake forward. (Illustration 34.)
- Pull up the hand brake to release the wheels. (Illustration 35.)



Illustration 34.



Illustration 35.

- In order to disassemble/assemble the wheels, push the wheel cover and pull the wheel.
(Illustration 36., 37.)



Illustration 36.



Illustration 37.

7.8. Replacement tires

Before replacing the tires, refer to the following notes.

First, let out the air in the innertube and then by pressing on the tire towards the center of the rim to loosen contact tire-rim. Using a special tool, pull the tire outside the rim, which will allow you to translate the whole tire to the outside. After this procedure, removing the innertube is very natural.

Before installing a new inner tube, please note that:

- the tape on the trousers was not broken and did not have any foreign objects,
- the exterior edge of the wheels was smooth, without foreign objects and other damage.

If necessary, the bottom tape and rim must be cleaned thoroughly.

In order to mount the tires, apply the bottom tape with the hole on the valve, and then put it in the rim. Keeping that position, should put the bottom tape around the perimeter of the rim, paying attention to whether there has been turning of the tape or whether it has been damaged, and well covered the tensioners and spokes.

Pulling the tire edge of the rim, insert the valve in the rim slot. Inflate slightly the innertube so that it obtained a round shape and then put the innertube in the rim. To be absolutely sure, check that the innertube is not twisted (in case of breakdowns or deform slightly inflate the air). By using your hand put the tire edge on the rim starting from the opposite side of the valve and ending at the valve.

Always check both sides of the wheel, do the innertube has not been pressed by the tire and the tires properly adheres to the rim. The valve must gently push inward, and then squeeze outside in order to well-settled in the slot.

To check the correctness mounting the tires, should inflate enough air to be able to press a thumb to the rim . Then check the tire adheres equally on both sides of the rim and whether it is centrally placed. When there defect should, let out the air from the wheel and improve innertube. When the innertube is placed properly, should inflate the tire to maximum pressure and screw the cap on the tip of the valve.



WARNING !

- when using the wrong tools or the wrong steps can lead to serious damage to wheels and tires;
- in case of application of the plastic rim does not apply the bottom tape ;
- pay attention in order to during the installation between the rim and the botton tape does not allow any foreign objects that may damage the component or hurt hand
- Observe proper level of tire inflation in accordance with the values provided by the manufacturer on the tire edge;
- should use a pump with a pressure gauge or manual pump.

7.9. Batery charging.


- Place the wheelchair close to a standard electric al outlet.
- Check if all the wires are connected properly.
- The colour of LED light of battery condition meter will change from yellow to green when the batteries are fully charged.



WARNING !

- Do not expose the wheelchair and battery to high temperature and sources of fire, especially while charging the batteries. Not following this rule may lead to ignition and explosion.
- Always charge the batteries in a well ventilated place.
- Never touch the charger , batteries with wet hands or other objects. Keep the charger, wires and batteries away from water, rain or damp conditions. Not following the above rules may cause electric shock.
- Always charge the batteries in min 0: and max 40: temperature. Not following this rule may lead to battery damage.
- Do not interrupt charging till the battery condition meter lights green.
- If the battery condition meter does not light, disconnect the wires for 10 min and try to connect them again. If it does not help, consult the dealer.
- If you charge the batteries for more than 12 hours but the battery condition meter is still yellow, please, contact the dealer.

Power supply parameters	Entrance 230v, 2a	exit Dc 24v, max. 5a 144w
----------------------------	-----------------------------	---------------------------------



7.10. Problem solving



Illustration 38.

- No voltage indication when power is on. Check the wires connection.
- Open the black cover placed at the back of the frame.
- Disconnect the plug and connect them to the socket to start charging.
- The LED light of control panel will flash during charging. The green colour indicates that charging has been completed and the batteries are charged, orange colour- batteries are charged in 50 %, red colour – batteries are charged in 30 %.

7.11. Batteries

- The wheelchair is equipped with deep-cycle batteries which are sealed and maintenance free.
- Please, remember that some situations require more energy (e.g. going up the incline, overloading, difficult weather conditions).
- In order to protect the batteries and make them work longer, please, follow the rules below:
 - Do not allow the batteries to discharge completely. Always charge the batteries after daily use.
 - While charging the batteries for the first time, disconnect the plug after 40 hours, when the batteries are fully charged.
 - If you do not use the wheelchair regularly, charge the batteries every three days.



WARNING !

- Never remove the caps from sealed batteries.
- The battery acid may cause blindness or burns. If the battery acid contact eyes, skin, clothes or any other object which has contact with skin, rinse with water as soon as possible and consult your doctor.
- Always use only Hero batteries –using other batteries may lead to wheelchair damage.
- Always make sure that the wheelchair and batteries are properly protected.
- Please, remember that the batteries will degenerate and the distance range will decrease after some time. If you notice that your batteries are less efficient, you should replace them with new one.
- Please, remember that for optimal effect, you should replace both batteries at the same time.



WARNING !

- Please, remember that incorrect connection may damage the batteries.
- If the batteries were expose to low temperatures (if the battery is frozen), do not try to charge them!
- Cold or frozen batteries should be kept for couple of days in room temperature in order to warm them up and charge.

7.12. Transport

In particular, attention should be paid to:

- Manual standing wheelchair should be transported without the user. Manual standing wheelchair is not intended as an additional seat in the transport vehicle.
- Straps should be fastened to the strengths and permanent elements of the vehicle;
- Manual standing wheelchair should be fasten only for the frame;
- After fastened, the manual standing wheelchair, should not move in any direction :
- During transport manual standing wheelchair the brakes should be activated;

7.13. Maintenance

If you follow the maintenance rules outlined in this section, you can help to ensure that your wheelchair will give you years of trouble-free operation. If you have any doubt as to your wheelchair care or operation, contact the local dealer.

ANNUAL MAINTENANCE

Annual maintenance will help you to check if the wheelchair is in good condition and if it works properly. It will also help to prevent any possible complications which may arise from any malfunctions.

The annual maintenance should include:

- Checking the tires.
- Checking the batteries and all connections.
- Checking if the VSI control panel works correctly.
- Checking the frame.
- Checking the upholstery.



WARNING !

- Protect the VSI control panel against shocks and blows.
- Do not expose the wheelchair to high/low temperatures and humidity.
- Keep the VSI control panel clean.
- Protect your wheelchair against contact with water.
- Keep the frame and upholstery clean. Clean them with soft, dry cloth, never use any abrasive. Do not use any detergents directly to VSI control panel, batteries or electrical parts.
- Grease the wheels every three months.
- Use only the batteries provided with the wheelchair. If you want to replace them, contact the dealer.

- Charge the batteries regularly. Always make sure that the charger is properly plugged into the socket. Never unplug the charger by pulling the wire.
- Check all the wires. If you notice any damage, secure it immediately.
- Store the wheelchair in a dry place, humidity may damage the frame and electronic al system.



Usage period - 5 years since production date

EVERYDAY USE AND MAINTANCE

- Clean the upholstery and frame with a moist cloth, without the addition of detergents and other similar chemicals.
- Systematically remove the dirt from the movable elements of the wheelchair (screws, resistance mechanisms).
- Keep clean all wheelchair elements for to save functional parameters of the product. Storing in rooms of very high humidity (baths, laundries etc.) or in the open air, in the combination with the lack of hygiene will cause the fast loss of functional and aesthetic parameters.

You must not use:

- solvents,
- toilet detergents,
- sharp brushes and hard objects,
- detergents that contain chlorine,
- corrosive detergents,
- Do not clean the exerciser with pressure or steam aggregates!


Repair

- In case of reparation, please contact the producer.








Customer service

- If you have any questions or need help, contact the producer.

8. Technical paramets

Model	Baby Hero
motor	6000 N
batteries	2.9Ah×12V×CPS
batteries weight	2.8 kg
charger	2A
brakes	Hamulec ręczny
max loading	100 kg 
Admissible tyre pressure	4,5 bar

Net Weight (depending on the model) BABY HERO: K, S,	30,5 kg
---	---------

Rozmiar								Przednie koło (cm)	Tylne koło (cm)
K	1045	555	860	475	340	340	370	15,24	55,80
S	1060	585	900	550	400	380	420	20,32	55,80

9. Symbols



Confirmation of compliance with EU standards



Manufacturer



Date of production



Read the user's manual



Warning




Catalog number



Serial number



Protect from humidity



Protect from sunlight



Permitted user weight.



Use inside and outside buildings



The product is intended to transport by car.



The product is not designed for use as a seat to transport in a motor vehicle.



The product is intended to the transport by plane.



The possibility of fingers jamming.



Do not stand on the footrest

10. Contact details



mdh sp. z o.o.

22/24 ks. W. Tymienieckiego street 90-349 Łódź, Poland

Tel. +48 42 674 83 84

Fax. +48 42 636 52 21

www.mdh.pl www.viteacare.com