

Primero Steel and aluminium wheelchair

Catalogue number: VCWK9C NEW

VITEA[®]
CARE

PRODUCT CHARACTERISTICS:

- foldable, aluminum frame of cross construction,
- adjustable angle of backrest position and seat depth
- detachable nylon black upholstery,
- Velcro® straps allow to adjust backrest depth,
- thick, soft cushion,
- safety belt,
- lifted, removable, angle adjustable footrests and calf strap,
- adaptable position of front and rear wheel axle,
- lifted, detachable armrests,
- adjustable position of the armrests in the horizontal and vertical plane
- adjustable handle height,
- rear barrier which provides additional stabilization and facilitates folding the wheelchair,
- rear anti-tip wheels of adjustable height,
- full front wheels, size 200x50,
- pneumatic quick-release rear wheels 24",
- pump,
- brakes of adjustable strength,
- aluminium handrims.



USE:

General use wheelchair, at home and for strolls. Suitable for patients after leg (legs) amputation, for patients during post-operative period.

	16"	18"
Seat width	40 cm	45 cm
Seat depth	40 – 52,5 cm	40 – 52,5 cm
Wheelchair width	61 cm	68 cm
Wheelchair height	91,5 – 103 cm	91,5 – 103 cm
Wheelchair length	107 cm	107 cm
Wheelchair weight	19,67 kg	19, 67 kg
Maximum load	115 kg	115 kg