

Mobile Parapodium

model PJMP with table, PJM without table

PRODUCT CHARACTERISTICS:

- stable structure allows to ambulate without any help,
- available in two vest variants: with and without lateral supports, the vest without lateral supports is designed for more able patients, who do not need chest support,
- supplied with adjustable “pants” securing patient from falling down,
- leg separator,
- base equipped with swivel wheels with clamps,
- additional equipment: wrist holders, additional cushions , table with regulated height and tilt angle.

USE:

Rehabilitation device for patients with balance disorders and/or partially able lower limbs.



VITEA
CARE

	100	125	150	180	200
Vest width*	21 - 27 cm	25 - 36 cm	29 - 43 cm	37 - 61,5 cm	37 - 61,5 cm
Vest depth*	16,5 - 21 cm	20 - 26 cm	23 - 34 cm	26 - 38 cm	29 - 41 cm
Parapodium width (mobile base)	36 - 46 cm	40 - 51 cm	53 - 69 cm	66 - 88 cm	66 - 88 cm
Parapodium length (mobile base)	76 cm	77 cm	80 cm	80 cm	80 cm
Parapodium height (mobile base)	64 - 82 cm	72 - 91 cm	89 - 116 cm	115 - 152 cm	125 - 158 cm
Height	85 - 105 cm	105 - 125 cm	125 - 155 cm	155 - 185 cm	185 - 205 cm
Maximum weight support	20 kg	35 kg	70 kg	90 kg	110 kg

* measurements without upholstery

