

Combo Mobile Parapodium

model PJMC

PRODUCT CHARACTERISTICS:

- consists of one vest and two bases: static and mobile,
- adjustable to patient's needs,
- allows the patient to be safe and comfortable in a vertical position,
- additional equipment: velcro straps for feet, reinforced knee mount connectors, wrist holders, wedge, additional cushions, table with regulated height and tilt angle.

USE:

Rehabilitation device for patients with balance disorders and/or partially able lower limbs.



VITEA[®]
CARE

	100	125	150	180	200
Vest width*	21 - 27 cm	25 - 36 cm	29 - 43 cm	37 - 61,5 cm	37 - 61,5 cm
Vest depth*	16,5 - 21 cm	20 - 26 cm	23 - 34 cm	26 - 38 cm	29 - 41 cm
Parapodium width (static base)	35 - 43 cm	40 - 53 cm	50 - 62 cm	59 - 83 cm	59 - 83 cm
Parapodium length (static base)	67 - 87 cm	67 - 87 cm	77 - 101 cm	77 - 101 cm	77 - 101 cm
Parapodium height (static base)	67 - 85 cm	80 - 99 cm	98 - 124 cm	123 - 160 cm	133 - 170 cm
Parapodium width (mobile base)	36 - 46 cm	44 - 51 cm	53 - 69 cm	66 - 88 cm	66 - 88 cm
Parapodium length (mobile base)	76 cm	77 cm	80 cm	80 cm	80 cm
Parapodium height (mobile base)	64 - 82 cm	72 - 91 cm	89 - 116 cm	115 - 152 cm	125 - 158 cm
Height	85 - 105 cm	105 - 125 cm	125 - 155 cm	155 - 185 cm	185 - 205 cm
Maximum weight support	20 kg	35 kg	70 kg	90 kg	110 kg

* measurements without upholstery

