

Anti-bedsore alternating pressure bubble mattress

catalogue number VCM202

PRODUCT CHARACTERISTICS:

- alternating pressure bubble mattress made of PCV,
- smooth manual pressure adjustment, adjusted individually to patient's weight,
- working cycle 12 minutes,
- silent electric pump,
- facings which can be fixed to an ordinary mattress,
- colour: beige.

USE:

Prevention, and treatment support of bedsores.

| | 202D | 202B |
|---------------------|--------|--------|
| Length | 200 cm | 200 cm |
| Width | 90 cm | 90 cm |
| Height | 6,5 cm | 7 cm |
| Max. weight support | 120 cm | 120 cm |



VITEA[®]
CARE

