

Bath bench

catalogue number DRVW042

PRODUCT CHARACTERISTICS:

- aluminium structure,
- seat is made of easy cleanable polyurethane,
- flow openings,
- irregular adjustment of handles width.

USE:

The bench facilitates patients transfer from and into the bath bench. It also facilitates every day toilet activities for people with motor limitations. It also provides a safe position for patient with bending limitations in hips e.g. after the hip endoprosthesis.

Adjustment of leg stops width 62 cm, 70 cm



VITEA[®]
CARE



mdh **CE**